



NAVINA FOR
CHILDREN

Securing freedom
for families

Navina™ Systems

Navina Smart:

Intelligent and precise

Navina™ Smart has been developed with specific focus on safety, efficiency and predictability. It comes with an electronic control unit and an app where sessions and progress can be monitored. The electronic control unit ensures that the therapy is consistent in its performance every time, no matter who performs the procedure.

Control, safety and predictability


- Electronic control unit (no manual pumping required)
- Personalized settings, consistent every time
- Built-in safety margins - maximum balloon size, water volume and water flow rate
- Automatic rejection of used catheters for safety and hygiene



Navina Smart app: The digital bowel diary

The Navina Smart app is a good tool to monitor treatment progress and satisfaction. The app stores data and settings from the Navina Smart unit after each completed irrigation, and pairs this data with the user's satisfaction rating. For you as a healthcare professional it means more control, efficient follow-up meetings, quality data for progress analysis and treatment recommendations.

- Stores treatment data and settings (balloon size, water amount, water flow rate, duration time)
- Allows the child to rate each irrigation creating involvement in the procedure
- Designed to help find optimal treatment parameters faster
- Intuitive and easy to use for everyone



Navina Classic:

The intuitive and convenient way to perform TAI manually

Navina™ Classic can be a complement to the Smart System, or it is simply for those who prefer manual pumping.

- Easy pumping for both air (balloon) and water
- Flexible and durable⁷ balloon
- Sturdy water container with double scale – for optimal visibility
- Intuitive – color-coding on both tubing and control unit (connectors and pumps)
- Symbols – easy to understand



Catheter or cone catheter

The Navina™ Systems works with both catheter and cone.

- Hydrophilic catheter suitable for children
- Flexible and durable⁷ balloon, seals the rectum throughout treatment
- Hydrophilic cone designed to fit snugly against the bottom
- Soft and flexible material with good grip



Securing freedom for families.

Bowel dysfunction is a condition that can have many different causes. Nonetheless the result is the same. A life turned upside down, time taken away from the more fun things in life, not only for the child but for the whole family.

Child:

- Constipation, fecal incontinence, and other bowel issues.
- Embarrassment, social insecurity.
- Limited time for play, activities and friends.


Caregivers:

- Stress and worries.
- Time-consuming care, limited time for rest of the family.

TAI is a long-term solution for children – and their families

Transanal irrigation (TAI) is a therapy with the potential to help many children living with bowel problems. It is a safe, effective and controlled way to empty the bowel regularly using only normal water.

- Effective treatment in children¹
- Reduces fecal incontinence and fecal soiling^{2,3,4}
- Improve family satisfaction with the child's bowel management¹
- Improves Quality of Life^{2,3,4,5}
- Reduces time spent on bowel management³
- Improves patient independence^{6,7}



1. Mosiello G et al. J Pediatr Gastroenterol Nutr. Mar 2017;64(3):343-352

References

1. Mosiello G et al. J Pediatr Gastroenterol Nutr. Mar 2017;64(3):343-352
2. Bray L and Sanders C. Spinal Cord. Feb 2013;51(2):88-93
3. Midrio P et al. Colorectal Dis. Jan 2016;18(1):86-93
4. Choi EK et al. Spinal Cord Nov 2014; 1-5
5. Ausili E et al. Spinal Cord. Jul 2010;48(7):560-565
6. Wide et al. Acta Paediatr. Nov 2014;103(11):1159-1164
7. Lopez Pereira P et al. J Pediatr Urol. Apr 2010;6(2):134-138

At Wellspect we value the people behind our success as a leading provider of life-changing products for bowel and bladder management such as the well-known brands LoFric® and Navina™. From the thousands of users and healthcare professionals worldwide who inspire our innovative solutions, we know that working together is the best way to advance continence care, giving our users more time for life. Building on over 30 years of life-improving performance, we passionately strive to make a difference every day to everyone who needs our products and services. Wellspect. A Real Difference.

wellspect.us

