



Intended use

LoFric® Hydro-Kit™ single use urinary catheter is intended for intermittent urinary catheterization.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

Packaging

30 catheters per customer box. 2 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

Plastic packaging: Polyethylene butylacrylate (EBA).

Product range

Reference	Type	Length	Size
421 08 XXX*	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 08
421 10 XXX*	LoFric® Hydro-Kit™ Nelaton (paediatric)	20 cm	CH 10
423 08 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 08
423 10 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 10
423 12 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 12
423 14 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 14
423 16 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 16
423 18 XXX*	LoFric® Hydro-Kit™ Nelaton (female)	20 cm	CH 18
420 08 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 08
420 10 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 10
420 12 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 12
420 14 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 14
420 16 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 16
420 18 XXX*	LoFric® Hydro-Kit™ Nelaton (male)	40 cm	CH 18
425 10 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 10
425 12 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 12
425 14 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 14
425 16 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 16
425 18 XXX*	LoFric® Hydro-Kit™ Tiemann (male)	40 cm	CH 18

*Generic article number without the three digit suffix specific for a region or country destination when distributing an article.

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

Wellspect. A real difference.

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

