



Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>1</sup>

Smooth eyelets.



Foldable into handy, discreet pocket size.



User friendly customer box.

### Intended use

LoFric® Origo™ single use urinary catheter is intended for intermittent urinary catheterization.

### Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

### Packaging

30 catheters per customer box. 4 customer boxes per shipping box.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sachet:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Item package:** Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Bottom foil: Polyolefin.

### Product range

Reference	Type	Length	Size
442 08 00	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 08
442 10 00	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 10
442 12 00	LoFric® Origo™ Nelaton (paediatric)	30 cm	CH 12
443 10 00	LoFric® Origo™ Nelaton (male)	40 cm	CH 10
443 12 00	LoFric® Origo™ Nelaton (male)	40 cm	CH 12
443 14 00	LoFric® Origo™ Nelaton (male)	40 cm	CH 14
443 16 00	LoFric® Origo™ Nelaton (male)	40 cm	CH 16
443 18 00	LoFric® Origo™ Nelaton (male)	40 cm	CH 18
444 10 00	LoFric® Origo™ Tiemann (male)	40 cm	CH 10
444 12 00	LoFric® Origo™ Tiemann (male)	40 cm	CH 12
444 14 00	LoFric® Origo™ Tiemann (male)	40 cm	CH 14
444 16 00	LoFric® Origo™ Tiemann (male)	40 cm	CH 16
444 18 00	LoFric® Origo™ Tiemann (male)	40 cm	CH 18

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

**Wellspect. A real difference.**

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)



Red Leaf Medical Inc.  
17-3250 Ridgeway Drive,  
Mississauga, ON L5L 5Y6,  
Canada  
phone:1-877-563-7422  
[www.redleafmedical.com](https://www.redleafmedical.com)

