



*Discreet design.*

*User-friendly handle for good grip and non-touch technique.*

*Adhesive tab on back for reliable hanging.*

*Instant activation - sterile salt solution included, squeeze and it's ready to use.*

*Peel-open tab.*

*The outer packaging doubles as a hygienic and discreet disposal pouch after use.*

*Not manufactured with PVC, phthalates or latex.*

*Smooth eyelets.*



*LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.<sup>1</sup>*

### Intended use

LoFric® Sense™ single use urinary catheter for females is intended for intermittent urinary catheterization.

### Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

### Packaging

30 catheters per customer box. 4 customer boxes per shipping box.

### Material composition

**Hydrophilic surface coating:** Polyvinyl pyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sachet:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Item package:** Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

### Product range

Reference	Type	Length	Size
416 08 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 08
416 10 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 10
416 12 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 12
416 14 XX*	LoFric® Sense™ Nelaton (female)	15 cm	CH 14

*\*Generic article number without the two digit suffix specific for a region or country destination when distributing an article.*

*The effective lengths is 10 cm. Effective length, length of the catheter that can be inserted into the body.*



*User friendly customer box.*

**LoFric® Sense™**

1. Waller et al 1997; Dunger et al. 1988

At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and healthcare professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

**Wellspect. A real difference.**

For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)

