

- Peel-open tab.
- Adhesive tab on back for reliable hanging.
- Instant activation - sterile salt solution included, squeeze and it's ready to use.
- Adjustable Insertion Grip allows better grip and control, promoting a non-touch technique.

Not manufactured with PVC, phthalates or latex.



LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{1,2}

• Smooth eyelets.^{3,4}



Foldable into handy, discreet pocket size.



Nelaton Coudé

Product description

LoFric® Origo™, single use, hydrophilic urinary catheter with salt solution for activation.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

Package

30 catheters per customer box. 4 customer boxes per shipping box.

Material composition

Hydrophilic surface coating: Polyvinyl pyrrolidone (PVP).

Catheter: Polyolefin-based elastomer (POBE).

Sachet: Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).


Item packaging:

Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).

Bottom foil: Polyolefin.

Product range

| Reference | Type | Length | Size |
|-----------|------------------------------------|--------|-------|
| 442 08 40 | LoFric® Origo™ Nelaton (pediatric) | 12" | 08 FR |
| 442 10 40 | LoFric® Origo™ Nelaton (pediatric) | 12" | 10 FR |
| 442 12 40 | LoFric® Origo™ Nelaton (pediatric) | 12" | 12 FR |
| 443 10 40 | LoFric® Origo™ Nelaton (male) | 16" | 10 FR |
| 443 12 40 | LoFric® Origo™ Nelaton (male) | 16" | 12 FR |
| 443 14 40 | LoFric® Origo™ Nelaton (male) | 16" | 14 FR |
| 443 16 40 | LoFric® Origo™ Nelaton (male) | 16" | 16 FR |
| 443 18 40 | LoFric® Origo™ Nelaton (male) | 16" | 18 FR |
| 444 10 40 | LoFric® Origo™ Coudé (male) | 16" | 10 FR |
| 444 12 40 | LoFric® Origo™ Coudé (male) | 16" | 12 FR |
| 444 14 40 | LoFric® Origo™ Coudé (male) | 16" | 14 FR |
| 444 16 40 | LoFric® Origo™ Coudé (male) | 16" | 16 FR |
| 444 18 40 | LoFric® Origo™ Coudé (male) | 16" | 18 FR |



User-friendly packaging.

LoFric® Origo™



References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.us](https://www.wellspect.us).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

[wellspect.us](https://www.wellspect.us)

Wellspect HealthCare, 1235 Friendship Rd, Suite 205, Braselton, GA 30517
Toll-free: 1 877 456 3742 . www.wellspect.us

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only