



Foldable into handy, discreet pocket size.



Flexible Coudé

## Product description

LoFric® Origo™, single use, hydrophilic urinary catheter with salt solution for activation.

## Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on package.

## Package

30 catheters per customer box.  
4 customer boxes per shipping box.

## Material composition

**Hydrophilic surface coating:** Polyvinylpyrrolidone (PVP).

**Catheter:** Polyolefin-based elastomer (POBE).

**Sachet:** Laminate foil of polyethylene terephthalate (PET), aluminum and polyethylene (PE). Sterile water with sodium chloride (NaCl).

**Sleeve:** Polyethylene (PE)

## Item packaging:

Top foil: Laminate foil of polyethylene terephthalate (PET) and polyethylene (PE).  
Bottom foil: Polyolefin.

## Product range

Reference	Type	Length	Size
449 12 40	LoFric® Origo™ Sleeve Coudé FX (male)	16"	FR 12
449 14 40	LoFric® Origo™ Sleeve Coudé FX (male)	16"	FR 14
449 16 40	LoFric® Origo™ Sleeve Coudé FX (male)	16"	FR 16



User-friendly packaging.



**LoFric® Origo™**

## References

1. Waller et al., Spinal Cord, 1997
2. Dunger et al., Acta Endocrinologica, 1988
3. Lundgren et al., Spinal Cord, 2000
4. Newman et al., Springer, 2018

At Wellspect, we develop innovative products and solutions that make a real difference in the lives of people managing continence care needs, empowering them to live each day to the fullest. We have been a trusted leader in healthcare for over 40 years with our product brands LoFric<sup>®</sup> and Navina<sup>™</sup>, and we continue to expand access to our user-friendly solutions. Every day, we strive to maintain the highest quality standards while minimizing the environmental impact of our products. By working closely with users and healthcare professionals, we aim to improve clinical outcomes and enhance quality of life for people around the world.

**Wellspect. A Real Difference.**

For more information about our products and services, please visit [Wellspect.us](https://Wellspect.us).

Join the conversation on LinkedIn, X, Facebook and Instagram.

**wellspect.us**

Wellspect HealthCare, 1235 Friendship Rd, Suite 205, Braselton, GA 30517.

Toll-free: +1 877 456 3742 .

**LoFric**



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.  
Phone: +46 31 376 40 00. [www.wellspect.com](https://www.wellspect.com)

**CE** 2797

**R<sub>x</sub>** Only