

Product description

LoFric® Elle™, a single use hydrophilic, ready to use urinary catheter, intended for intermittent urinary catheterization.

Sterilization and shelf life

Sterilized by irradiation. Use before expiry date on the package.

Package

30 catheters per customer box.
2 customer boxes per shipping box.
5 catheters per sample pack.

Material composition

Hydrophilic surface coating:
Polyvinyl pyrrolidone (PVP).

Catheter:
Polyolefin-based elastomer (POBE).

Wetting fluid:
Water-based wetting fluid with glycerol.

Item packaging:
Polypropylene.

Product range

Reference	Type	Length	Size
418 10 50	LoFric® Elle™ Nelaton (female)	4" *	10 FR
418 12 50	LoFric® Elle™ Nelaton (female)	4" *	12 FR
418 14 50	LoFric® Elle™ Nelaton (female)	4" *	14 FR
418 10 05	LoFric® Elle™ Nelaton (female)	4" *	10 FR <small>sample pack</small>
418 12 05	LoFric® Elle™ Nelaton (female)	4" *	12 FR <small>sample pack</small>
418 14 05	LoFric® Elle™ Nelaton (female)	4" *	14 FR <small>sample pack</small>


*Effective length, length of the catheter that can be inserted into the body.



Sample pack, 5 catheters.



User-friendly packaging.



- Designed to be ergonomic, user-friendly, with an angulated handle – for good grip and non-touch technique. Connects to a urine collection bag or extension tube.
- Ready to use – ready when you need it, to enable fast handling and catheterization.
- Slim and discreet design.
- Lower container is recyclable as plastic (polypropylene).
- Not manufactured with PVC, phthalates or latex.
- Smooth eyelets.^{1,2}
- LoFric catheters are coated with the unique Urotonic™ Surface Technology. This ensures the catheter remains smooth and safe, both during insertion and withdrawal.^{3,4}



LoFric® Elle™

References

1. Lundgren et al., Spinal Cord, 2000
2. Newman et al., Springer, 2018
3. Waller et al., Spinal Cord, 1997
4. Dunger et al., Acta Endocrinologica, 1988

At Wellspect we develop innovative continence care solutions that improve quality of life for people with bladder and bowel problems. We inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 40 years with our product brands LoFric® and Navina™. We always aim to minimize the environmental impact of our products and passionately strive to become climate neutral. We work together with users and healthcare professionals to improve clinical outcome in a sustainable way, now and for the future.

Wellspect. A Real Difference.

For more information about our products and services, please visit [Wellspect.us](https://www.wellspect.us).

Join the conversation on LinkedIn, Twitter, Facebook and Instagram.

wellspect.us

Wellspect HealthCare, 1235 Friendship Rd, Suite 205, Braselton, GA 30517
Toll-free: 1-877-456-3742 . www.wellspect.us

LoFric



Manufacturer: Wellspect HealthCare, Aminogatan 1, P.O. Box 14, SE-431 21 Mölndal, Sweden.
Phone: +46 31 376 40 00. www.wellspect.com

CE 2797

Rx Only